



## Human CTRP-1 Antibody

Mouse monoclonal

Catalog # M2003 (unconjugated)

M2003-b (biotin)

100 µg/ 200 µg (biotin)

Lot # L10003

---

Description	
Immunogen	Recombinant Human Complement C1q tumor necrosis factor-related protein 1 (CTRP1) (Met1-Pro281, with a C-terminal DDK and His tag expressed in HEK 293) Accession # NP_112230; Uniprot # Q9BXJ1
Reactivity	Human, not tested for other species
Clone	SLH-154
Purification	Protein A or G purified from hybridoma culture supernatant
Applications	ELISA
Formulation	Lyophilized from 1X PBS with Trehalose Reconstitute in 100µL ddH <sub>2</sub> O to 1 mg/mL (unconjugated) Reconstitute in 200µL ddH <sub>2</sub> O to 1 mg/mL (biotin)
Shipping	Ambient or cold
Storage	5 years at -20°C to -80°C 1 month at 4°C after reconstitution with preservative
Expiration	1/28/2021

For Research Use Only