


Description	
Source	HEK293 expressed Met1-Ser298, with a C-terminal DDK and His tag Accession # NP_001176
Purity	>95% by SDS-PAGE
N-terminal Sequence	Gln21, predicted
Predicted Molecular Mass	32.1kDa
SDS-PAGE	41-45 kDa, reducing conditions
Formulation	Lyophilized from 1X PBS with Trelahose
Preparation and Storage	
Reconstitution	Reconstitute at 50 µg/mL in sterile PBS containing at least 0.1% protein
Shipping	Ambient or cold
Storage	5 years at -80°C; 2 years at -20°C or 4°C Upon reconstitution: 5 years at -80°C; 2 years at -20°C
Expiration	1/28/2021



For Research Use Only