


Interleukin-2 (IL-2) is a cytokine signaling molecule that regulates white blood cells. It is produced primarily by T cells after activation by an antigen. IL-2 signals through the IL-2 receptor which leads to the activation of at least 2 major signaling pathways which mediate cell growth, survival, activation induced cell death and differentiation. A pro-inflammatory cytokine, IL-2 activates T cells to produce TNFα as well as IFNγ and controls immune response by contributing to the development of regulatory T cells. Therefore IL-2 plays an important role in antigen specific immune responses

Description	
Source	HEK293, Met1-Gln155, with C-terminal His Tag
UniProt	B3VBJ5
Purity	>99% by SDS-PAGE
Predicted Molecular Mass	17.6 kDa
SDS-PAGE	17.6 kDa
Formulation	Lyophilized from 1X PBS with Trehalose
Preparation and Storage	
Reconstitution	Reconstitute at 25 µg/mL in sterile PBS containing at least 0.1% protein
Shipping	Ambient
Storage	5 years at -80°C; 2 years at -20°C or 4°C Upon reconstitution: 5 years at -80°C; 2 years at -20°C
Expiration	See lot specific CoA
<div> <div> 118kDa- 98kDa- 62kDa- 49kDa- 38kDa- 28kDa- 14kDa- 6kDa- </div>  </div>	
Hamster IL-2 (2µg) on a Coomassie stained 4-12% SDS-PAGE gel, reducing conditions.	

For Research Use Only

AMPERSAND BIOSCIENCES, LLC

6459 State Route 30, Lake Clear, NY 12945 | PHONE: (518) 891 0204 | info@ampersandbio.com

Rev.1