


Tumor Necrosis Factor alpha (TNF α) is Potent cytokine involved in systemic inflammation. Abundantly released following trauma or infection it is vital in the production of pro-inflammatory cytokine cascade.

Description	
Immunogen	Recombinant Ferret TNF α expressed in E. coli. Val90-Leu234.
Reactivity	Ferret
Source	Mouse monoclonal IgG ₁
Purification	Protein G purified from hybridoma cell culture supernatant
Applications	ELISA (Detection)
Formulation	Lyophilized from PBS with Trehalose Reconstitute in 100 μ L ddH ₂ O to 1 mg/mL
Shipping	Ambient
Storage	2 years at -20°C to -80°C as supplied 1 month at 4°C after reconstitution with preservative 1 year at -20°C to -80°C after reconstitution
Expiration	2 years from date of receipt
 <p>Ferret TNFα Luminex Assay</p> <p>Capture: Ferret TNFα clone FTT-3H12-2G11 (Cat # M2049) Detection: Ferret TNFα clone FT17-2C12 (Cat # M2059) Antigen: Recombinant Ferret TNFα [E. coli]</p>	

For Research Use Only

Ampersand Biosciences products may not be sequenced, reverse engineered, or investigated in any way to determine composition without the expressed written consent of Ampersand Biosciences.