


**Interferon Gamma-Induced Protein 10 (IP-10, CXCL10)** is a chemokine that is secreted during an immune response, specifically following IFN gamma stimulation. Following activation by IFN gamma IP-10 initiates Th1-cells and if uncontrolled can lead to a cytokine storm and uncontrolled inflammation.

Description	
Immunogen	Full length recombinant ferret IP-10.
Reactivity	Ferret
Source	Mouse monoclonal IgG <sub>1</sub>
Purification	Protein G purified from hybridoma cell culture supernatant
Applications	ELISA (Capture)
Formulation	Lyophilized from PBS with Trehalose Reconstitute in 100µL ddH <sub>2</sub> O to 1 mg/mL
Shipping	Ambient
Storage	2 years at -20°C to -80°C as supplied 1 month at 4°C after reconstitution with preservative 1 year at -20°C to -80°C after reconstitution
Expiration	2 years from date of receipt
<div style="display: flex; justify-content: space-between;"> <div style="width: 45%;">  <p style="font-size: small;">             Ferret IP-10              SPL (Log): R<sup>2</sup> = 0.99815              Active (●) Inactive (◻) Curve fit (—)           </p> </div> <div style="width: 50%;"> <p><b>Ferret IP-10 Luminex Assay</b></p> <p>Capture: Ferret IP-10 clone <b>FTIP10-4H8</b> (Cat # M2048)            Detection: Ferret IP-10 clone FTIP10-4E7 (Cat # M2057)            Antigen: Recombinant Ferret IP-10</p> </div> </div>	

For Research Use Only

Ampersand Biosciences products may not be sequenced, reverse engineered, or investigated in any way to determine composition without the expressed written consent of Ampersand Biosciences.